

This certificate is granted and awarded by the authority of the Nadcap Management Council to:

Liebherr Aerospace Toulouse SAS

408 Avenue Des Etats Unis Bp 2010 Toulouse Cedex, 31016 France

This certificate demonstrates conformance and recognition of accreditation for specific services, as listed in www.eAuditNet.com on the Qualified Manufacturer's List (QML), to the revision in effect at the time of the audit for:

Materials Testing Laboratories

Certificate Number: 6642220518 Expiration Date: 28 February 2025 Accreditation Length: 18 Months

Jay Solomond

Executive Vice President & Chief Operating Officer



SCOPE OF ACCREDITATION

Materials Testing Laboratories

Liebherr Aerospace Toulouse SAS 408 Avenue Des Etats Unis Bp 2010 Toulouse Cedex, 31016 France

This certificate expiration is updated based on periodic audits. The current expiration date and scope of accreditation are listed at: www.eAuditNet.com - Online QML (Qualified Manufacturer Listing).

In recognition of the successful completion of the PRI evaluation process, accreditation is granted to this facility to perform the following:

AC7000 Rev A - AUDIT CRITERIA FOR NADCAP ACCREDITATION

AC7101/1 Rev G - Nadcap Audit Criteria for Materials Testing Laboratories – General Requirements for All Laboratories (to be used on audits BEFORE 10-Dec-2023)

AC7101/3 Rev D - Nadcap Audit Criteria for Materials Testing Laboratories – Mechanical Testing (to be used on audits on/after 4 December 2016)

- (A) Room Temperature Tensile
- (A1) Room Temperature Tensile with Elastic (Young's) Modulus
- (O) High Cycle Fatigue

Lab Type - Lab Type

Captive

t-frm-0004 29-May-20